
p64 MW Flow Modulation Device

phenox

INSTRUCTIONS FOR USE

2	DEUTSCH GEBRAUCHSANWEISUNG	29	DANSK BRUGSANVISNING	56	БЪЛГАРСКИ ИНСТРУКЦИИ ЗА УПОТРЕБА
5	ENGLISH INSTRUCTIONS FOR USE	32	NORSK BRUKSANVISNING	59	EESTI KASUTUSJUHEND
8	FRANÇAIS MODE D'EMPLOI	35	POLSKI INSTRUKCJA UŻYTKOWANIA	62	HRVATSKI UPUTE ZA UPORABU
11	ITALIANO ISTRUZIONI PER L'USO	38	SLOVENŠČINA NAVODILA ZA UPORABO	65	SLOVENČINA NÁVOD NA POUŽITIE
14	ESPAÑOL INSTRUCCIONES DE USO	41	ČEŠTINA NÁVOD K POUŽITÍ	68	LIETUVIŲ K. NAUDOJIMO INSTRUKCIJA
17	PORTUGUÊS INSTRUÇÕES DE UTILIZAÇÃO	44	ROMANĂ INSTRUCȚIUNI DE UTILIZARE	71	TÜRKÇE KULLANMA TALİMATLARI
20	NEDERLANDS GEBRUIKSAANWIJZING	47	MAGYAR HASZNÁLATI UTASÍTÁS	74	РУССКИЙ ЯЗЫК ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
23	SVENSKA BRUKSANVISNING	50	ΕΛΛΗΝΙΚΑ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	77	اللغة العربية إرشادات الاستخدام
26	SUOMI KÄYTTÖOHJEET	53	LATVISKI LIETOŠANAS PAMĀCĪBA		

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Hersteller

phenox GmbH
Lise-Meitner-Allee 31
44801 Bochum, Deutschland
Telefon: +49 234 36 919 0
Telefax: +49 234 36 919 19

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1 x **p64 MW** (HPC) Flow Modulation Device (in den Ausführungen **p64 MW** oder **p64 MW HPC**)

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- in whom antiplatelet and/or anticoagulation therapy is contraindicated or antiplatelet therapy did not start in a timely manner prior to treatment,
- in whom angiography demonstrates the anatomy is not appropriate for endovascular treatment, such as severe vessel tortuosity or stenosis.

Package content

1 x **p64 MW** (HPC) Flow Modulation Device (of version **p64 MW** or **p64 MW HPC**)

Product description


Fig. 1: **p64 MW** (HPC) implant and delivery system in introducer sheath


Fig. 2: Delivery system and detached **p64 MW** (HPC) implant

Please note that in the following text the term **p64 MW** (HPC) stands for both device versions, **p64 MW** (uncoated) and **p64 MW HPC** (coated).

The **p64 MW** (HPC) Flow Modulation Device is a tubular vascular implant that consists of 64 interwoven Nitinol wires

⑦ which are filled with a platinum core to ensure visibility under X-ray fluoroscopy.

The HPC coating (HPC: Hydrophilic Polymer Coating) of version **p64 MW HPC** covers the entire implant ⑦ and reduces the initial adherence of thrombocytes and lowers so the risk of thrombus generation.

The delivery system ④ has a platinum marker ⑥ at the distal end of the transport tube ④ and another one at the distal wire tip ⑧ to allow the operator to determine its position.

The attachment of the implant ⑦ to the delivery system ④ follows the friction locking principle: The proximal end of the implant

⑦ is secured between a soft polymer pad (at distal end of the transport tube ④) and an introducer sheath ⑤ (after removal of the sheath, the function is performed by the microcatheter), in such a way that pushing and pulling of the implant ⑦ is enabled.

The product is stored in an introducer sheath ⑤ and is transferred into a microcatheter with an inside diameter of 0.021 inches (0.53 mm). This sheath ⑤ is moved proximally during insertion of **p64 MW** (HPC) to enable complete passage through the microcatheter.

A white FluorSAFE Marker ⑩ on the transport tube ④ identifies the position to which the device can be advanced inside the microcatheter without the device tip leaving the microcatheter.

The implant ⑦ self-expands as it leaves the microcatheter. Until it has been fully deployed in the target vessel, the implant ⑦ can be completely recovered back into the microcatheter to allow it to be repositioned or removed. The point of maximum implant deployment that allows for implant recovery is indicated by a platinum marker ⑥ at the distal end of the transport tube ④: As long as the marker ⑥ is located inside the microcatheter the implant ⑦ can be completely recovered.

The **p64 MW** (HPC) is always deployed by means of a coordinated movement, whereby the microcatheter is withdrawn and the delivery system ④ is advanced to avoid any movement of the distal implant end by the shortening effect. Due to the shortening effect, the distal delivery wire tip ⑧ moves distally during deployment. To counteract this movement in order to avoid, e.g., the entry of the delivery wire tip ⑧ into distal sensitive vessels, the delivery wire tip ⑧ can be moved proximally after the torquer ② is released prior to the implant ⑦ being completely deployed. To do this the white torquer ② at the proximal end of the delivery system ④ is loosened and replaced by any standard torquer (compatible with a 0.014 inch or 0.016 inch (0.36 or 0.41 mm) microguidewire); this torquer is then locked more proximally to the end of the delivery wire ①. The delivery wire ① is then withdrawn out of the transport tube ④. The transport tube ④ has an additional handle ③ at its proximal end for easier handling.

The implant ⑦ is always recovered back into the microcatheter by means of a coordinated movement, whereby the microcatheter is advanced and the delivery system is withdrawn.

All manipulations are carried out under X-ray visualization. After final control of deployment and position, the implant ⑦ is completely deployed and detached from the delivery system ④ by withdrawal of the microcatheter.

Intended use

The **p64 MW** (HPC) Flow Modulation Device is a self-expanding, tubular vascular implant and allows the controlled and selective modulation of blood flow in extra- and intracranial arteries. In addition, the physical properties of the **p64 MW** (HPC) straighten the target vessel slightly and reinforce it.

These properties aid the endovascular reconstruction of diseased arteries along their cervical and intracranial course.

Indications

The **p64 MW** (HPC) Flow Modulation Device is a self-expanding, tubular implant and is used in the endovascular treatment of vascular diseases such as

- saccular and fusiform aneurysms and pseudoaneurysms,
- vascular dissections in the acute and chronic phases and
- vascular perforations and AV fistulae.

Contraindications

Treatment is contraindicated in patients

Compatibility

All **p64 MW** (HPC) models are compatible with the Rebar-18 microcatheter (Medtronic, USA) which has an inner diameter of 0.021 inches (0.53 mm).


In its relaxed state, the diameter of the **p64 MW** (HPC) is 0.4 mm larger than the nominal diameter. The length specifications on the packaging describe the clinically usable length.

The **p64 MW** (HPC) must be used according to the specifications regarding the minimum and maximum target vessel diameters which are stated on the packaging.

p64 MW (HPC) is available in the following versions:

- Uncoated: **p64 MW** (REF no. P64-MW-XXX-XX)
- Coated with HPC: **p64 MW HPC** (REF no. P64-MW-HPC-XXX-XX)

The size specifications are indicated by the REF no., and are also stated on the packaging:

Information on size selection

- Select diameter so that the deployed diameter comes as close as possible to the target vessel diameter, in order to achieve proper vessel wall apposition.
- Do not use the implant in target vessels whose diameter is not within the range of application specified on the packaging.
- Caution: Substantial oversizing (selection of a **p64 MW** (HPC) with a range of application considerably above the diameter of the target vessel) poses the risk of incorrect deployment (incomplete expansion).
- Caution: Undersizing (selection of a **p64 MW** (HPC) with a range of application below the diameter of the target vessel) leads to insufficient fixation of the **p64 MW** (HPC) within the vessel and allows blood to flow around the outside of the implant (a so-called "endoleak"). The implant is then unstable, is subject to migration and hemodynamically ineffective.
- Ensure that the implant overlaps the lesion distally and proximally. If the selected product is too short or too long, it can be removed and replaced with a suitable one.
- Ensure that the implant does not end proximally in a narrow vessel curve because this may constrain a full proximal expansion. Choose an implant length which results in a complete coverage of the proximal vessel curve by the **p64 MW** (HPC).

Information on the selection of patients and lesions

If compliance to the antiplatelet medication cannot be guaranteed following the implantation of a **p64 MW** (HPC), thrombotic closure of the implant and the vessel around it may occur within just a few days. Patients who cannot comply with the prescribed medication may not be suitable for treatment with a **p64 MW** (HPC).

From the time of implantation of a **p64 MW** (HPC), several weeks or months may pass before an aneurysm is no longer a risk. In this period, no complete protection from a (fresh) rupture/bleed can be guaranteed. Therefore, patients who are in the acute phase post aneurysm rupture should be treated by options that offer greater protection from re-rupture/bleed.

Medication

Prior and following implantation of a **p64 MW** (HPC) antiplatelet medication is necessary as described in the chapter "Recommended procedure". Be mindful of possible interactions with other medications (e.g. with proton pump inhibitors, Ibuprofen, Metamizole).

In vitro test results and anecdotal clinical experience demonstrate that the version **p64 MW HPC** can provide a reduced surface thrombogenicity. In justified exceptional cases the reduced thrombogenicity can allow the implantation under single antiplatelet medication, only if there no reasonable alternative therapy given. Here particular attention is to be paid to at least three days medication prior to treatment. The achieved platelet inhibition is more intensive by using P2Y12 inhibitors (Prasugrel, Ticagrelor) than by using ASA. For safety reasons, the efficacy of the antiplatelet medication is always to be verified by means of appropriate tests (e.g. Multiplate, VerifyNow, PFA).

Single antiplatelet medication may have an increased risk of thromboembolic events if multiple devices have been implanted in a telescoping fashion. The risk of thrombus formation may be increased after subarachnoid hemorrhage, after trauma, during pregnancy, after major surgery, during inflammatory diseases, fever, thrombocytosis.

In general dual antiplatelet medication in the context of flow diverter implantation is safer than mono medication concerning the risk of thromboembolic events. Dual antiplatelet medication, however, carries a higher risk of hemorrhagic events.

ASA is less efficacious than P2Y12 inhibitors concerning the protection from thromboembolic events. Several conditions increase the required ASA dosage significantly (intracranial hemorrhage, pregnancy, trauma, surgery, thrombocytosis, fever, pneumonia...). The action of ASA is antagonized by Ibuprofen and Metamizole. ASA is in several countries available as a variant which can be given IV. ASA does usually not cause hemorrhagic issues if surgery is required.

Prasugrel has been reported to prevent thrombus formation on HPC coated devices. These are so far anecdotal observations. Controlled trials are pending. The risk of hemorrhagic complications from Prasugrel remains a concern.

Ticagrelor might be a compromise for single antiplatelet medication. The short-acting time requires consistent intake.

Recommended procedure

Preparation of procedure and patient, platelet aggregation inhibition and patient testing

1. Gather and document as complete a case history as possible, especially regarding the history of the present illness, comorbidities, previous interventions and current medication.
2. As far as possible, inform the patient and document the patient's consent to the planned intervention, pointing out the possible complications and potential consequences (disability, care dependency or death). In cases where patients cannot give consent themselves, their relatives should, as far as

possible, be asked if they know what the presumed patient's will would be. Otherwise, in case of emergency, the rules of emergency care for incapacitated patients apply, subject to different institutional or national requirements.

- Promptly take all measures necessary for suitable pretreatment with drugs to ensure the inhibition of platelet aggregation.**

Based on current knowledge, dual antiplatelet therapy when implanting the **p64 MW** (HPC) and similar products is suitable for the prevention of thrombus formation caused by the implant. For this purpose, 100 mg of ASA and 75 mg of Clopidogrel can be administered orally every day for at least 3 days before the planned intervention. Alternatively, 500 mg of ASA and 600 mg of Clopidogrel can be administered orally as one-off doses on the day before treatment.

Preloading with high dosage might be less reliable than loading with the regular dosage for several days in terms of protection against thrombus formation. High-dosage preloading may result in hyper-response which may cause hemorrhagic complications (e.g., intracerebral and subarachnoid hemorrhage).

In vitro test results and anecdotal clinical experience demonstrate that the version **p64 MW** HPC can provide a reduced surface thrombogenicity. If required and justified by individual circumstances, the reduced thrombogenicity of **p64 MW** HPC may allow the implantation under single antiplatelet medication. It is recommended to discuss this procedure with the patient and his or her legal representatives. In this case, particular attention has to be paid to a pretreatment for at least three days prior to the treatment. The achieved platelet inhibition is stronger by using P2Y12 inhibitors (Prasugrel, Ticagrelor) than by using ASA.

If ASA is used as a single medication, 2 x 100 mg ASA PO daily (1-0-1) is recommended. If a P2Y12 receptor inhibitor is used, Prasugrel is apparently more efficacious than Clopidogrel and Ticagrelor. Prasugrel might increase the risk of hemorrhagic complications if compared to Clopidogrel. If Ticagrelor is used the short-acting time of this drug has to be considered. The regular dosage is 2 x 90 mg Ticagrelor PO daily (1-0-1). Inconsistent intake of Ticagrelor is associated with an increased risk of thromboembolic events.

The safety of the treatment is increased if the effective inhibition of platelet function is verified by means of an appropriate test (e.g. Multiplate, VerifyNow, PFA) before the intervention. Regarding substitutes in case of resistance to Clopidogrel and the use of Gp IIb/IIIa antagonists, we refer you to the respective current scientific publications. See also chapter "Medication".
- A CT or MRI examination of the cranium and, where necessary, the throat is advisable beforehand, in order to ensure a comprehensive preliminary diagnosis.
- The diagnostic angiography and endovascular treatment should be carried out under general anesthesia with neuromuscular relaxation and invasive hemodynamic monitoring. During anesthesia, aim to maintain suitable systolic blood pressure values.
- After preparing both groins, a 6F or 8F catheter is inserted, preferably into the right femoral artery.
- Then moderate **heparinization** should begin, which should also last for the duration of the intervention. An intravenous dose of between 3000 and 5000 units of heparin has proven suitable in practice. Where available, determining the ACT ("activated clotting time") is advisable.
- Angiographic visualization of the internal and external carotid arteries on both sides and of the vertebral artery on at least one side is recommended, along with the respective dependent vessels. Enlarged images and, where necessary, oblique images of the affected vessel(s) are recommended.
- The target vessel(s) for the endovascular treatment must be defined.
- A 6F guide catheter or the combination of an 8F guide catheter and a suitable extension catheter or distal access catheter is inserted into the afferent cervical vessel, taking steps to avoid vasospasm.
- It is important that the implant is only inserted into target vessels of suitable size.**

Measure the diameter of the target vessel, where the distal and proximal ends of the **p64 MW** (HPC) are to be anchored, as accurately as possible.

Carefully observe and respect the specifications regarding the minimum and maximum vessel diameters on the packaging, as well as the instructions regarding selection of a model of the correct size (see Information on size selection). The length of the **p64 MW** (HPC) must be selected so that the implant covers the lesion at the distal and proximal ends by at least a few millimetres.

Introduction of microcatheter

- Never probe against resistance!**

Insert a suitable microcatheter with a corresponding microguidewire into the target vessel using a hemostatic valve and pressurized irrigation. Here, the use of so-called "road map" technology is advisable. Aim to position the tip of the microcatheter 10–15 mm distal to the treatment target. Once the target vessel treatment segment has been reached, carefully pull on the microcatheter in order to remove any excess catheter length and straighten the catheter.
 - Remove the microguidewire from the microcatheter under X-ray fluoroscopy.
- #### Preparation and introduction of p64 MW (HPC)
- Take the sterile device in its dispenser snail out of the packaging. Release the proximal end of **p64 MW** (HPC) and pull it together with introducer sheath out of the dispenser snail.
 - With the aid of a tight-closing hemostatic valve and under continuous pressurized irrigation with heparinized physiological saline solution, the **p64 MW** (HPC) is transferred from its introducer sheath into the microcatheter. For this purpose, the hemostatic valve is opened. The introducer sheath of the **p64 MW** (HPC) is inserted through the open valve. The hemostatic valve is closed carefully and the introducer sheath of the **p64 MW** (HPC) is flushed by retrograde entry of the irrigation fluid.
 - Once the introducer sheath of the **p64 MW** (HPC) is completely flushed in this manner, it is advanced carefully until it reaches the distal end of the hub adapter of the microcatheter. The introducer sheath is held fixed in this position. The **p64 MW** (HPC) is then advanced from the introducer sheath into the microcatheter using the transport tube to which the implant is fixed. This process is continued until the white Fluorosafe Marker of the transport tube reaches the proximal end of the introducer sheath.
 - The introducer sheath is then pulled back proximally as far as the handle on the transport tube. (During the further procedure the sheath stays on the delivery system.)

The **p64 MW** (HPC) is pushed further forward until the Fluorosafe Marker of the Transport tube reaches the entry of the hemostatic valve. This procedure does not need fluoroscopy because the Fluorosafe Marker identifies the position to which the device can be advanced without the device tip leaving the microcatheter.

The process of inserting the **p64 MW** (HPC) generally corresponds to that of inserting other similar implants. Should you encounter particular resistance which can only be overcome with effort, the implant and possibly also the microcatheter must be removed and the vessel accessed once more.
 - Never push the p64 MW (HPC) delivery system tip beyond the distal tip of the microcatheter. This can lead to a dissection or perforation of the target vessel.**

The **p64 MW** (HPC) is slowly advanced to the tip of the microcatheter under continuous fluoroscopy. The distal tip of the delivery system should reach the tip of the microcatheter.

Deployment of p64 MW (HPC)

- Release the implant fully by carefully and very slowly withdrawing the microcatheter up to point that the implant can still be recovered back into the microcatheter. The point of maximum implant deployment that allows for implant recovery is indicated by a platinum marker at the distal end of the transport tube: As long as marker is located inside the microcatheter the implant can be completely recovered. Once the distal end of the implant is fully expanded and anchored in the distal vessel, continue to deploy the implant by continuously pushing on the delivery system in order to facilitate the expansion of the **p64 MW** (HPC). In order to ensure optimal wall apposition, the deployment must be a coordinated effort of continuous pushing of the delivery system and adjustments (advancing or withdrawing) of the microcatheter so that the microcatheter is centralized longitudinally along the vessel. The release of the **p64 MW** (HPC) should take place under fluoroscopy in order to ensure that the implant is properly deployed and the distal end has not moved.

Repositioning of distal delivery wire tip (optional)

- Note that the distal delivery wire tip moves distally during implant deployment!**

To counteract this movement and to avoid, e.g., the entry of the delivery wire tip into distal sensitive vessels, the delivery wire tip can be moved proximally after the torquer is released while the implant is not completely deployed. To do this the white torquer at the proximal end of the delivery system is loosened and replaced by any standard torquer (compatible with a 0.014 inch or 0.016 inch (0.36 or 0.41 mm) microguidewire); this torquer is then locked more proximally to the end of the delivery wire. The delivery wire is then withdrawn out of the transport tube. The transport tube has an additional handle at its proximal end for easier handling.

Continuation of deployment

- The **p64 MW** (HPC) is self-expanding and, when properly deployed, apposes itself against the vascular wall. The implant may over-expand at the neck of the aneurysm due to the increased diameter at that point. Correct deployment can be verified by visualizing the platinum-filled braiding wires of the implant.
- Injecting approx. 6–10 ml of X-ray contrast medium through the guide catheter allows one to check whether the aneurysm/dissection/target vessel has been satisfactorily covered by the deployment and release of the **p64 MW** (HPC).
- If the radial deployment of the **p64 MW** (HPC) is insufficient or the position or the model size selected is unsuitable, the implant can be recovered into the microcatheter, if the distal marker of the transport tube is still inside the microcatheter, in order to allow the implant to be repositioned, redeployed or completely removed. If the delivery wire tip was moved proximally before, it must be ensured that the distal wire tip is placed again distally to the distal compressed implant end and the white torquer is locked again on the transport tube.

For repositioning or removal the microcatheter is advanced while the delivery system is slowly withdrawn.

Detachment of p64 MW (HPC)

- Due to the radial expansion of the proximal end, a slight shortening of the implant takes place!**

If the position and deployment of the **p64 MW** (HPC) are satisfactory, the implant is **immediately** fully deployed and detached by complete withdrawal of the microcatheter.

The proximal implant end is thus exposed and it can fully expand.

When using DSA systems with a digital detector and CT technology ("flat panel detector CT", e.g. DynaCT [Siemens], XperCT, VasoCT [Philips]), the implant can be visualized on the sectional image. This has proven particularly effective in the evaluation of the deployment and apposition to the vessel wall.
- Remove the delivery system by gently withdrawing.
- Insufficient deployment of the **p64 MW** (HPC) can be improved by means of a subsequent balloon dilatation. As far as possible, the **p64 MW** (HPC) should appose against the vascular wall.

Implantation of another p64 MW (HPC)

- After the first **p64 MW** (HPC) is detached, if a subsequent telescoping device is required, gently advance the microcatheter through the **p64 MW** (HPC). When the microcatheter tip is distal to the **p64 MW** (HPC), gently retract the wire tip into the microcatheter and remove the delivery system completely out of the microcatheter. The microcatheter is now in position for a subsequent **p64 MW** (HPC) to be advanced and deployed.
- Injecting approx. 6–10 ml of X-ray contrast medium through the guide catheter allows one to check once more, if necessary, whether the target vessel has been sufficiently covered by the application of the **p64 MW** (HPC). This check should be repeated 10 to 15 minutes later where necessary.

Platelet aggregation inhibition and responder testing

- Take steps to ensure adequate inhibition of platelet aggregation. Proven medications following implantation include a 1 x 100 mg oral dose of ASA every day on an ongoing basis and a 75 mg oral dose of Clopidogrel every day for at least 12 months, but longer when necessary or on an ongoing basis. Be mindful of possible interactions with other medications (e.g. with proton pump inhibitors, Ibuprofen, Metamizole).
- In vitro test results and anecdotal clinical experience demonstrate that the version **p64 MW** HPC can provide a reduced surface thrombogenicity. In justified exceptional cases, the reduced thrombogenicity can allow the implantation under single antiplatelet medication, only if no reasonable alternative therapy is given. Here particular attention is to be paid to at least three days medication prior to treatment. The achieved platelet inhibition is more intensive by using P2Y12 inhibitors (Prasugrel, Ticagrelor) than by using ASA.
- For safety reasons, the efficacy of the antiplatelet medication is always to be verified by means of appropriate tests (e.g. Multiplate, VerifyNow, PFA).
- Single antiplatelet medication may have an increased risk of thromboembolic events if multiple devices have been implanted in a telescoping fashion. The risk of thrombus formation may be increased after subarachnoid hemorrhage, after trauma, during pregnancy, after major surgery, during inflammatory diseases, fever, thrombocytosis.
- See also chapter "Medication".

Precautions

- Microcatheters with other inner diameters (ID) than 0.021 inch (e.g., 0.017 or 0.027 inch) do not work at all. **p64 MW** (HPC) used in microcatheters with too large IDs lead to premature detachment of the implant inside the microcatheter.
- The **p64 MW** (HPC) may be deployed up to three (3) times in the target vessel. It must be considered that each deployment only occurs up to the point that the distal marker of the transport tube is still inside the microcatheter!

- For flushing place the introducer sheath of the **p64 MW** (HPC) inside the hemostatic valve of the microcatheter and flush it by the help of the connected irrigation fluid. Thorough flushing of the introducer sheath is essential in order to remove any trapped air bubbles.
- All manipulations must be carried out under fluoroscopic visualization.
- If the **p64 MW** (HPC) system is advanced beyond the distal end of the microcatheter, the vessel may be dissected or perforated.
- If the **p64 MW** (HPC) system can be advanced into the microcatheter only with great effort or navigated through the microcatheter only with great effort, remove the entire **p64 MW** (HPC) system out of the microcatheter as a precaution.
- Do not pull the deployed implant back through the vessel into the microcatheter. Instead, push the microcatheter over the **p64 MW** (HPC) while simultaneously fixating the delivery system to reposition and redeploy the implant if necessary.
- The **p64 MW** (HPC) is a delicate implant and requires careful handling. Never push the microcatheter onto the **p64 MW** (HPC) against resistance. Never twist the delivery system. If necessary, remove the **p64 MW** (HPC) along with the microcatheter.
- Forceful pulling or pushing or twisting on the delivery system may inadvertently detach the **p64 MW** (HPC) from the delivery system. In such a case, recovery with a foreign body retrieval device (e.g. Microsnare Kit, ev3/Covidien/Medtronic) is recommended.
- Extreme stress to the delivery system tip could lead to separation of some of its parts. In such cases, recovery with a foreign body retrieval device (e.g. Microsnare Kit, ev3/Covidien/Medtronic) is recommended.
- If there is the risk that the patient is unresponsive or not fully responding to dual antiplatelet therapy, the time between implant deployment and detachment is supposed to be as short as possible to avoid any insufficient deployment of proximal implant end by bonding of braiding wires caused by blood ingredients (e.g. fibrin).
Any additional treatment (e.g. coiling of the aneurysm while the associated microcatheter is "jailed" by the deployed **p64 MW** (HPC)) is supposed to be done after **p64 MW** (HPC) detachment.
- Do not detach the implant in case of poor wall apposition or if the implant is in a stretched condition. Improve expansion by re-deploying or removing/replacing the device.
- If not all the implant's braiding wires emerge from the soft pad at the distal end of the transport tube, careful movement of the delivery system helps to loosen them.
- Time between start of implant deployment and complete deployment resulting in a detachment must be as short as possible, in order to prevent any bonding effects by blood ingredients and ultimately poor proximal implant deployment.
- Insufficient deployment of the proximal implant end can be improved by manipulation with the delivery system, the microcatheter, or means of a subsequent balloon dilatation. If the proximal end of the **p64 MW** (HPC) does not open, recovery with a foreign body retrieval device (e.g. Microsnare Kit, ev3/Covidien/Medtronic) is recommended.
- If vasospasm is suspected in the affected vascular region, all necessary measures, e.g. medication, should be used to aid regression prior to implantation.
- For delivery wire withdrawal use a separate torque device locked to the delivery wire. A suitable torquer is compatible with microguidewires that are 0.014 or 0.016 inches (0.36 or 0.41 mm) in diameter.
- If the torquer cannot be removed from the transport tube (by unscrewing and pulling it in the proximal direction), unscrew the torquer device completely. A higher resistance may be felt prior to complete separation of torquer handle and cap. Remove the torquer handle while the torquer cap and inlay stay in place. Then continue using the separate torque device for the process of delivery wire withdrawal.
- If the implant cannot be detached or is not supposed to be detached after torquer release after the delivery wire was moved proximally and the **p64 MW** (HPC) implant needs to be removed, gently place the delivery wire again distally to enable that then the microcatheter can recover implant and delivery wire tip.
- If the distal marker of the transport tube displaces by the movement of the delivery wire, remove the entire **p64 MW** (HPC) system considering the above mentioned precautions.
- Patients with known hypersensitivity to nickel-titanium materials may suffer an allergic reaction to the implant.
- Certification: The implant may only be used by specialized and appropriately trained physicians. Completion of a product training course run by phenox GmbH is a prerequisite for the use of the **p64 MW** (HPC). At least three (3) interventions using the **p64 MW** (HPC) must be supervised by a physician or other qualified person delegated by phenox GmbH, and their course and outcome must be documented.

General information

- Keep away from heat. Store in a cool, dry place.
- Use is only permissible prior to the expiration date, as sterility is not guaranteed otherwise.
- Do not use if packaging is damaged, as sterility cannot be assumed otherwise.
- The device must be checked for damage before use. Do not use deformed or damaged devices, as function cannot be assumed otherwise.
- The product is intended for single use only. The implant must not be withdrawn into its introducer sheath, resterilized or reprocessed for use in other patients since it cannot be cleaned reliably.
- The delivery system and, where necessary, packaging components must be disposed of appropriately in marked containers.

Complications

The following complications, amongst others, may arise during or after implantation:

- Air embolism, embolism in distal vessels, vessel occlusion, thrombosis and cerebral ischemia
- Perforation, rupture, dissection and other arterial lesions
- Occlusion of side branch/perforator
- (Transient) stenosis of target vessel
- Vasospasm, occurrence of a pseudoaneurysm, intracranial hemorrhage
- Recurrence of aneurysm, retreatment of aneurysm
- Allergic reaction, infection



- Space-occupying infarction, neurological deficit including the consequences of a stroke
- Persistent vegetative state, death

Magnetic resonance imaging

Non-clinical tests have shown that the **p64 MW** (HPC) is suitable for MRI at a magnetic flux density of 3T. Under clinical conditions, 1.5T has proven to be unproblematic for the implant.

Symbols and their meaning

Symbols on the label:

	Caution
	Consult instructions for use
	Do not reuse
	Do not resterilize
	Batch code
	Sterilized using ethylene oxide
	Do not use if package is damaged
	Non-pyrogenic
	Use-by date
	Catalog number
	Content
	Keep dry. Keep away from sunlight.
	This product was brought onto the market in accordance with Directive 93/42/EEC concerning medical devices.
	Manufacturer
	Material: NiTi (Nickel Titanium alloy)
	Latex-free
	Phthalate-free

Limitation of liability

phenox GmbH shall not be liable for damage caused by use other than that intended of the product. phenox and p64 are registered trademarks of phenox GmbH in the Federal Republic of Germany and other states.

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phenox GmbH | Lise-Meitner-Allee 31 | D-44801 Bochum, Germany