Sažetak o sigurnosnoj i kliničkoj učinkovitosti (SSCP)

za skupinu medicinskih proizvoda klase III

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Tablica 10: Poznate publikacije o proizvodu za modulaciju protoka p64

Citat - p64 publikacije - status rujan 2024.
Publikacije samo o proizvodu p64 Aguilar Perez, M., et al., <i>Endovascular Treatment of Anterior Circulation Aneurysms With the p64 Flow Modulation Device: Mid- and Long-Term Results in 617 Aneurysms From a Single Center</i>. Oper Neurosurg (Hagerstown), 2021. 20(4): p. 355-363. Sirakov, S., et al., <i>The p64 Flow Diverter-Mid-term and Long-term Results from a Single Center</i>. Clin Neuroradiol, 2020. 30(3): p. 471-480. Hellstern, V., et al., <i>Microsurgical clipping and endovascular flow diversion of ruptured anterior circulation blood blister-like aneurysms</i>. Interv Neuroradiol, 2018. 24(6): p. 615-623. Aguilar Perez, M., et al., <i>In-stent Stenosis after p64 Flow Diverter Treatment</i>. Clin Neuroradiol, 2018. 28(4): p. 563-568. Morais, R., et al., <i>Endovascular treatment of intracranial aneurysms with the p64 flow diverter stent: mid-term results in 35 patients with 41 intracranial aneurysms</i>. Neuroradiology, 2017. 59(3): p. 263-269. Briganti, F., et al., <i>Mid-term and long-term follow-up of intracranial aneurysms treated by the p64 Flow Modulation Device: a multicenter experience</i>. J Neurointerv Surg, 2017. 9(1): p. 70-76. Briganti, F., et al., <i>p64 Flow Modulation Device in the treatment of intracranial aneurysms: initial experience and technical aspects</i>. J Neurointerv Surg, 2016. 8(2): p. 173-80. Fischer, S., et al., <i>Initial Experience with p64: A Novel Mechanically Detachable Flow Diverter for the Treatment of Intracranial Saccular Sidewall Aneurysms</i>. AJNR Am J Neuroradiol, 2015. 36(11): p. 2082-9. Hellstern, V., et al., <i>Endovascular Treatment of Posterior Circulation Saccular Aneurysms With the p64 Flow Modulation Device: Mid-and Long-Term Results in 54 Aneurysms From a Single Center</i>. Front Neurol, 2021. 12: p. 711863. De Beule, T., et al., <i>p64 flow diverter: Results in 108 patients from a single center</i>. Interv Neuroradiol, 2021. 27(1): p. 51-59. Bonafe, A., et al., <i>Diversion-p64: results from an international, prospective, multicenter, single-arm post-market study to assess the safety and effectiveness of the p64 flow modulation device</i>. J Neurointerv Surg, 2022. 14(9): p. 898-903. Publikacije o proizvodu p64 i drugim phenox preusmjerivačima protoka Vivanco-Suarez, J., et al., <i>Safety and efficacy of the p48 MW and p64 flow modulation devices: a systematic review and meta-analysis</i>. Neurosurg Focus, 2023. 54(5): p. E7. Bilgin, C., et al., <i>Phenox HPC and Phenox flow modulation devices for the endovascular treatment of intracranial aneurysms: a systematic review and meta-analysis</i>. J Neurointerv Surg, 2024. 16(7): p. 706-714. Hellstern, V., et al., <i>Flow diversion for unruptured MCA bifurcation aneurysms: comparison of p64 classic, p64 MW HPC, and p48 MW HPC flow diverter stents</i>. Frontiers in Neurology, 2024. 15. Publikacija o proizvodu p64 i sličnim preusmjerivačima protoka Cagnazzo, F., et al., <i>Treatment of Unruptured Distal Anterior Circulation Aneurysms with Flow-Diverter Stents: A Meta-Analysis</i>. AJNR Am J Neuroradiol, 2019. 40(4): p. 687-693. Zhou, G., et al., <i>Complications associated with the use of flow-diverting devices for cerebral aneurysms: a systematic review and meta-analysis</i>. Neurosurg Focus, 2017. 42(6): p. E17. Gory, B., et al., <i>Flow Diverters for Intracranial Aneurysms: The DIVERSION National Prospective Cohort Study</i>. Stroke, 2019. 50(12): p. 3471-3480. Bhogal, P., et al., <i>Treatment of Unruptured, Tandem Aneurysms of the ICA with a Single Flow Diverter</i>. Clin Neuroradiol, 2019. 29(4): p. 725-731. Wendl, C.M., et al., <i>Direct carotid cavernous sinus fistulae: vessel reconstruction using flow-diverting implants</i>. Clin Neuroradiol, 2017. 27(4): p. 493-501. Briganti, F., et al., <i>Postprocedural, midterm, and long-term results of cerebral aneurysms treated with flow-diverter devices: 7-year experience at a single center</i>. Neurosurg Focus, 2017. 42(6): p. E3.

Citat - p64 publikacije - status rujan 2024.

- Maybaum, J., et al., *Flow Diversion for Reconstruction of Intradural Vertebral Artery Dissecting Aneurysms Causing Subarachnoid Hemorrhage-A Retrospective Study From Four Neurovascular Centers*. Front Neurol, 2021. 12: p. 700164.
- Narata, A.P., et al., *Reversible Brain Edema Associated with Flow Diverter Stent Procedures: A Retrospective Single-Center Study to Evaluate Frequency, Clinical Evolution, and Possible Mechanism*. World Neurosurg, 2019. 122: p. e569-e576.
- Bhogal, P., et al., *Treatment of Unruptured, Saccular, Anterior Choroidal Artery Aneurysms with Flow Diversion: A Single Centre Experience*. Clin Neuroradiol, 2019. 29(3): p. 459-465.
- Yaltirik Bilgin, E., et al., *Endovascular Treatment of Intracranial Anterior Circulation Aneurysms with Flow Diverters: A Single Centre Experience with mid and long-term results*. Turk Neurosurg, 2017.
- Peschillo, S., et al., *Endovascular Treatment of Large and Giant Carotid Aneurysms with Flow-Diverter Stents Alone or in Combination with Coils: A Multicenter Experience and Long-Term Follow-up*. Oper Neurosurg (Hagerstown), 2017. 13(4): p. 492-502.
- Bhogal, P., et al., *The Use of Flow Diverting Stents to Treat Para-Ophthalmic Aneurysms*. Front Neurol, 2017. 8: p. 381.
- Bhogal, P., et al., *The Fate of Side Branches Covered by Flow Diverters-Results from 140 Patients*. World Neurosurg, 2017. 103: p. 789-798.
- Bhogal, P., et al., *Flow Diversion for the Treatment of MCA Bifurcation Aneurysms-A Single Centre Experience*. Front Neurol, 2017. 8: p. 20.
- Bhogal, P., et al., *Treatment of posterior circulation non-saccular aneurysms with flow diverters: a single-center experience and review of 56 patients*. J Neurointerv Surg, 2017. 9(5): p. 471-481.
- Guzzardi, G., et al., *Long-term follow-up in the endovascular treatment of intracranial aneurysms with flow-diverter stents: update of a single-centre experience*. Radiol Med, 2018. 123(6): p. 449-455.
- Bhogal, P., et al., *Management of Unruptured Saccular Aneurysms of the M1 Segment with Flow Diversion: A Single Centre Experience*. Clin Neuroradiol, 2018. 28(2): p. 209-216.
- Giorgianni, A., et al., *Flow Diversion for Acutely Ruptured Intracranial Aneurysms Treatment: A Retrospective Study and Literature Review*. J Stroke Cerebrovasc Dis, 2022. 31(3): p. 106284.
- Simgen, A., et al., *Endovascular treatment of unruptured intracranial aneurysms with flow diverters: A retrospective long-term single center analysis*. Neuroradiol J, 2023. 36(1): p. 76-85.
- Khanafar, A., et al., *Endovascular treatment of distal anterior cerebral artery aneurysms using flow modulation devices: mid- and long-term results from a two-center study*. Front Neurol, 2024. 15: p. 1368612.
- Abdel-Tawab, M., et al., *Efficacy and safety of flow diverters in posterior circulation aneurysms and comparison with their efficacy in anterior circulation aneurysms: A systematic review and meta-analysis*. Interv Neuroradiol, 2021. 27(5): p. 609-621.
- Alwakeal, A., et al., *Flow Diversion of Posterior Circulation Aneurysms: Systematic Review of Disaggregated Individual Patient Data*. AJNR Am J Neuroradiol, 2021. 42(10): p. 1827-1833.

Tablica 11: Poznate publikacije o proizvodu za modulaciju protoka p48 MW (HPC)

Citat - p48 MW (HPC) publikacije - status rujan 2024.

Publikacije samo o proizvodu p48 MW

- AlMatter, M., et al., *The p48 MW flow modulation device for treatment of unruptured, saccular intracranial aneurysms: a single center experience from 77 consecutive aneurysms*. CVIR Endovasc, 2020. 3(1): p. 39.
- Bhogal, P., et al., *The p48MW Flow Diverter-Initial Human Experience*. Clin Neuroradiol, 2021. 31(1): p. 135-145.

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- Vivanco-Suarez, J., et al., *Safety and efficacy of the p48 MW and p64 flow modulation devices: a systematic review and meta-analysis*. Neurosurg Focus, 2023. 54(5): p. E7.
- Bilgin, C., et al., *Phenox HPC and Phenox flow modulation devices for the endovascular treatment of intracranial aneurysms: a systematic review and meta-analysis*. J Neurointerv Surg, 2024. 16(7): p. 706-714.
- den Bergh, F.V., et al., *The p48 flow diverter: First clinical results in 25 aneurysms in three centers*. Interv Neuroradiol, 2021. 27(3): p. 339-345.
- Schob, S., et al., *Single-Center Experience With the Bare p48MW Low-Profile Flow Diverter and Its Hydrophilically Covered Version for Treatment of Bifurcation Aneurysms in Distal Segments of the Anterior and Posterior Circulation*. Front Neurol, 2020. 11: p. 1050.

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- Dabhi, N., et al., *Flow Diverter Devices for Treatment of Intracranial Aneurysms in Small Parent Vessels-A Systematic Review of Literature*. World Neurosurg, 2022. 162: p. 183-194.e7.
- Giorgianni, A., et al., *Flow Diversion for Acutely Ruptured Intracranial Aneurysms Treatment: A Retrospective Study and Literature Review*. J Stroke Cerebrovasc Dis, 2022. 31(3): p. 106284.
- Khanafer, A., et al., *Endovascular treatment of distal anterior cerebral artery aneurysms using flow modulation devices: mid- and long-term results from a two-center study*. Front Neurol, 2024. 15: p. 1368612.

Publikacije samo o proizvodu p48 MW HPC

- de Castro-Afonso, L.H., et al., *Treatment of distal unruptured intracranial aneurysms using a surface-modified flow diverter under prasugrel monotherapy: a pilot safety trial*. J Neurointerv Surg, 2021. 13(7): p. 647-651.
- de Castro-Afonso, L.H., et al., *Aspirin monotherapy in the treatment of distal intracranial aneurysms with a surface modified flow diverter: a pilot study*. J Neurointerv Surg, 2021. 13(4): p. 336-341.
- Bhogal, P., et al., *The p48_HPC antithrombogenic flow diverter: initial human experience using single antiplatelet therapy*. J Int Med Res, 2020. 48(1): p. 300060519879580.
- Aguilar-Perez, M., et al., *The p48 Flow Modulation Device with Hydrophilic Polymer Coating (HPC) for the Treatment of Acutely Ruptured Aneurysms: Early Clinical Experience Using Single Antiplatelet Therapy*. Cardiovasc Intervent Radiol, 2020. 43(5): p. 740-748.
- Pierot, L., et al., *Surface-modified flow diverter p48-MW-HPC: Preliminary clinical experience in 28 patients treated in two centers*. J Neuroradiol, 2021. 48(3): p. 195-199.
- de Castro-Afonso, L.H., et al., *Two year follow-up of distal unruptured intracranial aneurysms treated with a surface modified flow diverter under prasugrel monotherapy*. J Neurointerv Surg, 2023.

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- Vivanco-Suarez, J., et al., *Safety and efficacy of the p48 MW and p64 flow modulation devices: a systematic review and meta-analysis*. Neurosurg Focus, 2023. 54(5): p. E7.
- Bilgin, C., et al., *Phenox HPC and Phenox flow modulation devices for the endovascular treatment of intracranial aneurysms: a systematic review and meta-analysis*. J Neurointerv Surg, 2024. 16(7): p. 706-714.
- Lobsien, D., et al., *Aneurysm Treatment in Acute SAH with Hydrophilic-Coated Flow Diverters under Single-Antiplatelet Therapy: A 3-Center Experience*. AJNR Am J Neuroradiol, 2021. 42(3): p. 508-515.
- Guzzardi, G., et al., *Flow diverter stents with hydrophilic polymer coating for the treatment of acutely ruptured aneurysms using single antiplatelet therapy: Preliminary experience*. Interv Neuroradiol, 2020. 26(5): p. 525-531.
- Bhogal, P., et al., *Early clinical experience with the p48MW HPC and p64MW HPC flow diverters in the anterior circulation aneurysm using single anti-platelet treatment*. Interv Neuroradiol, 2022. 28(3): p. 266-276.
- Khanafer, A., et al., *Flow diversion with hydrophilic polymer coating with prasugrel as single antiplatelet therapy in the treatment of acutely ruptured intracranial aneurysms: a multicenter case series, complication and occlusion rates*. J Neurointerv Surg, 2024.
- den Bergh, F.V., et al., *The p48 flow diverter: First clinical results in 25 aneurysms in three centers*. Interv Neuroradiol, 2021. 27(3): p. 339-345.
- Schob, S., et al., *Single-Center Experience With the Bare p48MW Low-Profile Flow Diverter and Its Hydrophilically Covered Version for Treatment of Bifurcation Aneurysms in Distal Segments of the Anterior and Posterior Circulation*. Front Neurol, 2020. 11: p. 1050.
- Hellstern, V., et al., *Flow diversion for unruptured MCA bifurcation aneurysms: comparison of p64 classic, p64 MW HPC, and p48 MW HPC flow diverter stents*. Frontiers in Neurology, 2024. 15.

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- Ma, L., et al., *Flow Diverters with Surface Modification in Patients with Intracranial Aneurysms: A Systematic Review and Meta-Analysis*. World Neurosurg, 2024. 185: p. 320-326.e17.
- Schüngel, M.S., et al., *Distal Flow Diversion with Anti-Thrombotically Coated and Bare Metal Low-Profile Flow Diverters-A Comparison*. J Clin Med, 2023. 12(7).
- Gawlitza, M., et al., *A Systematic Literature Review and Meta-Analysis of the Treatment of Ruptured Intracranial Aneurysms with Hydrophilic Polymer and Phosphorylcholine-Coated Flow Diverters Under Single Antiplatelet Therapy*. World Neurosurg, 2023. 170: p. e791-e800.

Citat - p48 MW (HPC) publikacije - status rujan 2024.

- Monteiro, A., et al., *Treatment of ruptured intracranial aneurysms using the novel generation of flow-diverters with surface modification: A systematic review and meta-analysis*. Interv Neuroradiol, 2024. 30(3): p. 350-360.
- Maybaum, J., et al., *Flow Diversion for Reconstruction of Intradural Vertebral Artery Dissecting Aneurysms Causing Subarachnoid Hemorrhage-A Retrospective Study From Four Neurovascular Centers*. Front Neurol, 2021. 12: p. 700164.
- Schungel, M.S., et al., *Endovascular Treatment of Intracranial Aneurysms in Small Peripheral Vessel Segments- Efficacy and Intermediate Follow-Up Results of Flow Diversion With the Silk Vista Baby Low-Profile Flow Diverter*. Front Neurol, 2021. 12: p. 671915.
- Senol, Y.C., et al., *The safety profile of single antiplatelet therapy with flow diverters: Systematic review and meta-analysis*. Interv Neuroradiol, 2023: p. 15910199231168669.
- Goertz, L., et al., *Safety and efficacy of coated flow diverters in the treatment of ruptured intracranial aneurysms: a retrospective multicenter study*. J Neurointerv Surg, 2024.
- Khanafar, A., et al., *Endovascular treatment of distal anterior cerebral artery aneurysms using flow modulation devices: mid- and long-term results from a two-center study*. Front Neurol, 2024. 15: p. 1368612.

Tablica 12: Poznate publikacije o proizvodu za modulaciju protoka p64 MW (HPC)

Citat- p64 MW (HPC) publikacije - status rujan 2024.

Publikacije samo o proizvodu p64 MW HPC

- inters, H., et al., *First Experience of Three Neurovascular Centers With the p64MW-HPC, a Low-Profile Flow Diverter Designed for Proximal Cerebral Vessels With Antithrombotic Coating*. Front Neurol, 2021. 12: p. 724705.
- Petrov, A., et al., *Initial experience with the novel p64MW HPC flow diverter from a cohort study in unruptured anterior circulation aneurysms under dual antiplatelet medication*. Interv Neuroradiol, 2021. 27(1): p. 42-50.
- Hellstern, V., et al., *Use of a p64 MW Flow Diverter with Hydrophilic Polymer Coating (HPC) and Prasugrel Single Antiplatelet Therapy for the Treatment of Unruptured Anterior Circulation Aneurysms: Safety Data and Short-term Occlusion Rates*. Cardiovasc Intervent Radiol, 2022. 45(9): p. 1364-1374.
- Ernst, M., et al., *Multicenter study of the safety and effectiveness of intracranial aneurysm treatment with the p64MW-HPC flow modulation device*. Interv Neuroradiol, 2023: p. 15910199231220964.

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- Vivanco-Suarez, J., et al., *Safety and efficacy of the p48 MW and p64 flow modulation devices: a systematic review and meta-analysis*. Neurosurg Focus, 2023. 54(5): p. E7.
- Bilgin, C., et al., *Phenox HPC and Phenox flow modulation devices for the endovascular treatment of intracranial aneurysms: a systematic review and meta-analysis*. J Neurointerv Surg, 2024. 16(7): p. 706-714.
- Lobsien, D., et al., *Aneurysm Treatment in Acute SAH with Hydrophilic-Coated Flow Diverters under Single-Antiplatelet Therapy: A 3-Center Experience*. AJNR Am J Neuroradiol, 2021. 42(3): p. 508-515.
- Guzzardi, G., et al., *Flow diverter stents with hydrophilic polymer coating for the treatment of acutely ruptured aneurysms using single antiplatelet therapy: Preliminary experience*. Interv Neuroradiol, 2020. 26(5): p. 525-531.
- Bhogal, P., et al., *Early clinical experience with the p48MW HPC and p64MW HPC flow diverters in the anterior circulation aneurysm using single anti-platelet treatment*. Interv Neuroradiol, 2022. 28(3): p. 266-276.
- Khanafar, A., et al., *Flow diversion with hydrophilic polymer coating with prasugrel as single antiplatelet therapy in the treatment of acutely ruptured intracranial aneurysms: a multicenter case series, complication and occlusion rates*. J Neurointerv Surg, 2024.
- Hellstern, V., et al., *Flow diversion for unruptured MCA bifurcation aneurysms: comparison of p64 classic, p64 MW HPC, and p48 MW HPC flow diverter stents*. Frontiers in Neurology, 2024. 15.

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- Ma, L., et al., *Flow Diverters with Surface Modification in Patients with Intracranial Aneurysms: A Systematic Review and Meta-Analysis*. World Neurosurg, 2024. 185: p. 320-326.e17.
- Gawlitza, M., et al., *A Systematic Literature Review and Meta-Analysis of the Treatment of Ruptured Intracranial Aneurysms with Hydrophilic Polymer and Phosphorylcholine-Coated Flow Diverters Under Single Antiplatelet Therapy*. World Neurosurg, 2023. 170: p. e791-e800.

Citat- p64 MW (HPC) publikacije - status rujan 2024.

- Monteiro, A., et al., *Treatment of ruptured intracranial aneurysms using the novel generation of flow-diverters with surface modification: A systematic review and meta-analysis*. Interv Neuroradiol, 2024. 30(3): p. 350-360.
- Senol, Y.C., et al., *The safety profile of single antiplatelet therapy with flow diverters: Systematic review and meta-analysis*. Interv Neuroradiol, 2023: p. 15910199231168669.
- Goertz, L., et al., *Safety and efficacy of coated flow diverters in the treatment of ruptured intracranial aneurysms: a retrospective multicenter study*. J Neurointerv Surg, 2024.
- Khanafer, A., et al., *Endovascular treatment of distal anterior cerebral artery aneurysms using flow modulation devices: mid- and long-term results from a two-center study*. Front Neurol, 2024. 15: p. 1368612.

Bibliografija

1. Bhogal, P., et al., *Thrombogenicity of the p48 and anti-thrombogenic p48 hydrophilic polymer coating low-profile flow diverters in an in vitro human thrombin generation model*. Interv Neuroradiol, 2020. **26**(4): p. 488-493.
2. Lenz-Habijan, T., et al., *Hydrophilic Stent Coating Inhibits Platelet Adhesion on Stent Surfaces: Initial Results In Vitro*. Cardiovasc Intervent Radiol, 2018. **41**(11): p. 1779-1785.
3. Lenz-Habijan, T., et al., *Comparison of the Thrombogenicity of a Bare and Antithrombogenic Coated Flow Diverter in an In Vitro Flow Model*. Cardiovasc Intervent Radiol, 2020. **43**(1): p. 140-146.
4. Bannewitz, C., et al., *Evaluation of Antithrombogenic pHPC on CoCr Substrates for Biomedical Applications*. Coatings, 2021. **11**(1).
5. Fischer, S., et al., *Initial Experience with p64: A Novel Mechanically Detachable Flow Diverter for the Treatment of Intracranial Saccular Sidewall Aneurysms*. AJNR Am J Neuroradiol, 2015. **36**(11): p. 2082-9.
6. Aguilar Perez, M., et al., *Endovascular Treatment of Anterior Circulation Aneurysms With the p64 Flow Modulation Device: Mid- and Long-Term Results in 617 Aneurysms From a Single Center*. Oper Neurosurg (Hagerstown), 2021. **20**(4): p. 355-363.
7. Vivanco-Suarez, J., et al., *Safety and efficacy of the p48 MW and p64 flow modulation devices: a systematic review and meta-analysis*. Neurosurg Focus, 2023. **54**(5): p. E7.
8. Hellstern, V., et al., *Flow diversion for unruptured MCA bifurcation aneurysms: comparison of p64 classic, p64 MW HPC, and p48 MW HPC flow diverter stents*. Frontiers in Neurology, 2024. **15**.
9. AlMatter, M., et al., *The p48 MW flow modulation device for treatment of unruptured, saccular intracranial aneurysms: a single center experience from 77 consecutive aneurysms*. CVIR Endovasc, 2020. **3**(1): p. 39.
10. Ernst, M., et al., *Multicenter study of the safety and effectiveness of intracranial aneurysm treatment with the p64MW-HPC flow modulation device*. Interv Neuroradiol, 2023: p. 15910199231220964.
11. den Bergh, F.V., et al., *The p48 flow diverter: First clinical results in 25 aneurysms in three centers*. Interv Neuroradiol, 2021. **27**(3): p. 339-345.
12. De Beule, T., et al., *p64 flow diverter: Results in 108 patients from a single center*. Interv Neuroradiol, 2021. **27**(1): p. 51-59.
13. Schüngel, M.S., et al., *Distal Flow Diversion with Anti-Thrombotically Coated and Bare Metal Low-Profile Flow Diverters-A Comparison*. J Clin Med, 2023. **12**(7).
14. Khanafer, A., et al., *Flow diversion with hydrophilic polymer coating with prasugrel as single antiplatelet therapy in the treatment of acutely ruptured intracranial aneurysms: a multicenter case series, complication and occlusion rates*. J Neurointerv Surg, 2024.
15. Bonafe, A., et al., *Diversion-p64: results from an international, prospective, multicenter, single-arm post-market study to assess the safety and effectiveness of the p64 flow modulation device*. J Neurointerv Surg, 2022. **14**(9): p. 898-903.
16. Hellstern, V., et al., *Endovascular Treatment of Posterior Circulation Saccular Aneurysms With the p64 Flow Modulation Device: Mid-and Long-Term Results in 54 Aneurysms From a Single Center*. Front Neurol, 2021. **12**: p. 711863.
17. Sirakov, S., et al., *The p64 Flow Diverter-Mid-term and Long-term Results from a Single Center*. Clin Neuroradiol, 2020. **30**(3): p. 471-480.
18. Ma, L., et al., *Flow Diverters with Surface Modification in Patients with Intracranial Aneurysms: A Systematic Review and Meta-Analysis*. World Neurosurg, 2024. **185**: p. 320-326.e17.
19. Hellstern, V., et al., *Use of a p64 MW Flow Diverter with Hydrophilic Polymer Coating (HPC) and Prasugrel Single Antiplatelet Therapy for the Treatment of Unruptured Anterior Circulation Aneurysms: Safety Data and Short-term Occlusion Rates*. Cardiovasc Intervent Radiol, 2022. **45**(9): p. 1364-1374.
20. Winters, H., et al., *First Experience of Three Neurovascular Centers With the p64MW-HPC, a Low-Profile Flow Diverter Designed for Proximal Cerebral Vessels With Antithrombotic Coating*. Front Neurol, 2021. **12**: p. 724705.
21. Bhogal, P., et al., *The p48MW Flow Diverter-Initial Human Experience*. Clin Neuroradiol, 2021. **31**(1): p. 135-145.
22. Petrov, A., et al., *Initial experience with the novel p64MW HPC flow diverter from a cohort study in unruptured anterior circulation aneurysms under dual antiplatelet medication*. Interv Neuroradiol, 2021. **27**(1): p. 42-50.
23. Guzzardi, G., et al., *Flow diverter stents with hydrophilic polymer coating for the treatment of acutely ruptured aneurysms using single antiplatelet therapy: Preliminary experience*. Interv Neuroradiol, 2020. **26**(5): p. 525-531.

24. Bonafe, A., et al., *Diversion-p64: results from an international, prospective, multicenter, single-arm post-market study to assess the safety and effectiveness of the p64 flow modulation device*. J Neurointerv Surg, 2021.
25. Bilgin, C., et al., *Phenox HPC and Phenox flow modulation devices for the endovascular treatment of intracranial aneurysms: a systematic review and meta-analysis*. J Neurointerv Surg, 2023.
26. de Castro-Afonso, L.H., et al., *Two year follow-up of distal unruptured intracranial aneurysms treated with a surface modified flow diverter under prasugrel monotherapy*. J Neurointerv Surg, 2023.
27. King, R.M., et al., *Acute Thrombus Burden on Coated Flow Diverters Assessed by High Frequency Optical Coherence Tomography*. Cardiovasc Intervent Radiol, 2020. **43**(8): p. 1218-1223.
28. Shehata, M.A., et al., *Long-term outcomes of flow diversion for unruptured intracranial aneurysms: a systematic review and meta-analysis*. J Neurointerv Surg, 2023. **15**(9): p. 898-902.
29. Yarahmadi, P., et al., *Intra-procedural complications, success rate, and need for retreatment of endovascular treatments in anterior communicating artery aneurysms: a systematic review and meta-analysis*. Neurosurg Rev, 2022. **45**(5): p. 3157-3170.
30. de Castro-Afonso, L.H., et al., *Treatment of distal unruptured intracranial aneurysms using a surface-modified flow diverter under prasugrel monotherapy: a pilot safety trial*. J Neurointerv Surg, 2021. **13**(7): p. 647-651.
31. Bhogal, P., et al., *Early clinical experience with the p48MW HPC and p64MW HPC flow diverters in the anterior circulation aneurysm using single anti-platelet treatment*. Interv Neuroradiol, 2022. **28**(3): p. 266-276.
32. Luecking, H., et al., *Two- to five-year follow-up of 78 patients after treatment with the Flow Redirection Endoluminal Device*. Interv Neuroradiol, 2020. **26**(1): p. 38-44.
33. Kuhn, A.L., et al., *Impact of age on cerebral aneurysm occlusion after flow diversion*. J Clin Neurosci, 2019. **65**: p. 23-27.
34. Bhogal, P., et al., *Treatment of Unruptured, Saccular, Anterior Choroidal Artery Aneurysms with Flow Diversion : A Single Centre Experience*. Clin Neuroradiol, 2019. **29**(3): p. 459-465.
35. Schungel, M.S., et al., *Endovascular Treatment of Intracranial Aneurysms in Small Peripheral Vessel Segments-Efficacy and Intermediate Follow-Up Results of Flow Diversion With the Silk Vista Baby Low-Profile Flow Diverter*. Front Neurol, 2021. **12**: p. 671915.
36. Maus, V., et al., *The Use of Flow Diverter in Ruptured, Dissecting Intracranial Aneurysms of the Posterior Circulation*. World Neurosurg, 2018. **111**: p. e424-e433.
37. de Castro-Afonso, L.H., et al., *Aspirin monotherapy in the treatment of distal intracranial aneurysms with a surface modified flow diverter: a pilot study*. J Neurointerv Surg, 2021. **13**(4): p. 336-341.
38. Etminan, N., et al., *European Stroke Organisation (ESO) guidelines on management of unruptured intracranial aneurysms*. Eur Stroke J, 2022. **7**(3): p. V.
39. *Subarachnoid haemorrhage caused by a ruptured aneurysm: diagnosis and management [P] Evidence review for non-culprit aneurysms NICE guideline NG228 Subarachnoid haemorrhage Contents*. London: National Institute for Health and Care Excellence (NICE), 2022. NICE, 2022.
40. Zhao, J., et al., *Current Treatment Strategies for Intracranial Aneurysms: An Overview*. Angiology, 2018. **69**(1): p. 17-30.
41. Thompson, B.G., et al., *Guidelines for the Management of Patients With Unruptured Intracranial Aneurysms: A Guideline for Healthcare Professionals From the American Heart Association/American Stroke Association*. Stroke, 2015. **46**(8): p. 2368-400.
42. Henkes, H. and W. Weber, *The Past, Present and Future of Endovascular Aneurysm Treatment*. Clin Neuroradiol, 2015. **25 Suppl 2**: p. 317-24.
43. Lee, K.S., et al., *The evolution of intracranial aneurysm treatment techniques and future directions*. Neurosurg Rev, 2022. **45**(1): p. 1-25.
44. Papadopoulos, F., C.N. Antonopoulos, and G. Geroulakos, *Stent-Assisted Coiling of Unruptured Intracranial Aneurysms with Wide Neck*. Asian J Neurosurg, 2020. **15**(4): p. 821-827.
45. Phan, K., et al., *Meta-analysis of stent-assisted coiling versus coiling-only for the treatment of intracranial aneurysms*. J Clin Neurosci, 2016. **31**: p. 15-22.
46. Arimura, K. and K. Iihara, *Surgical Management of Intracranial Artery Dissection*. Neurol Med Chir (Tokyo), 2016. **56**(9): p. 517-23.
47. Kleindorfer, D.O., et al., *2021 Guideline for the Prevention of Stroke in Patients With Stroke and Transient Ischemic Attack: A Guideline From the American Heart Association/American Stroke Association*. Stroke, 2021. **52**(7): p. e364-e467.

48. Asif, K.S., et al., *Endovascular reconstruction for progressively worsening carotid artery dissection*. Journal of NeuroInterventional Surgery, 2015. **7**(1): p. 32-39.
49. Halbach, V.V., *Management of Vascular Perforations That Occur During Neurointerventional Procedures*. AJNR, 1991.
50. Poppert, *Zerebrale-Gefäßmalformationen*.
51. Steinmetz, H., *Leitlinien für Diagnostik und Therapie in der Neurologie. Kapitel Vaskuläre Erkrankungen-Subarachnoidalblutung (SAB), Stand 2012*. Deutsche Gesellschaft für Neurologie 2012.
52. Csippa, B., et al., *Hydrodynamic Resistance of Intracranial Flow-Diverter Stents: Measurement Description and Data Evaluation*. Cardiovasc Eng Technol, 2020. **11**(1): p. 1-13.
53. Brisman, J.L., J.K. Song, and D.W. Newell, *Cerebral Aneurysms*. The new england journal of medicine, 2006: p. 355:928-39.
54. Dholakia, R., et al., *Hemodynamics of Flow Diverters*. J Biomech Eng, 2017. **139**(2).
55. García-Núñez, J., et al., *A Study on the Efficacy and Safety of Pipeline Shield Stents and Nonmodified Surface Flow-Diverting Stents in Patients with Intracranial Aneurysms*. World Neurosurg, 2024. **183**: p. e781-e786.
56. Lieber, B.B., et al., *Particle image velocimetry assessment of stent design influence on intra-aneurysmal flow*. Ann Biomed Eng, 2002. **30**(6): p. 768-77.
57. Lieber, B.B., A.P. Stancampiano, and A.K. Wakhloo, *Alteration of hemodynamics in aneurysm models by stenting: influence of stent porosity*. Ann Biomed Eng, 1997. **25**(3): p. 460-9.
58. Kallmes, D.F., et al., *A new endoluminal, flow-disrupting device for treatment of saccular aneurysms*. Stroke, 2007. **38**(8): p. 2346-52.
59. Becske, T., et al., *Pipeline for uncoilable or failed aneurysms: results from a multicenter clinical trial*. Radiology, 2013. **267**(3): p. 858-68.